

# PL-Touch



## Introducing the Newest Line of Cold Lasers from Erchonia®

The PL-Touch laser from Erchonia features an easy to use graphical user interface (GUI) and four patented line generated lasers housed within two feather-weight removable energy delivery probes. With user definable protocols and a quick program mode, this handheld technology is easy to use and incredibly intuitive to operate. Experience the freedom of this incredible new technology— ***Touch your way to patient recovery!***

### The PL-Touch Advantage:

- *Touchscreen GUI*
- *Ergonomic design*
- *Made in the USA*

To learn more call or visit us online today!  
888-242-0571 | [www.erchonia.com](http://www.erchonia.com)

### SPECIFICATIONS

Configuration: (4) 7.5mW Line Generated Laser Diodes with Patented Optics

Wavelength: 635nm

Modulation: Constant Wave (CW) & Pulsed (Hz)

Display: Full Color TFT Touchscreen Technology for Ease of Use

Power Source: (2) 3.7 VDC Lithium-ion Polymer Batteries

Enclosure: Micro-Textured Anodized Aircraft Aluminum for Strength and Ease of Cleaning

Weight: Less than 2.0 lbs. (.91kg)

Accessories: Battery Charger and Carrying Case

Compliant to ISO 13485 Medical Device Quality, IEC 60601-1 Safety, 60601-1-2 EMC, IEC 60825-1 Laser Safety, CB Certified, CE Mark, CMDCAS

Laser Class 2 / Device Class II (USA); 2a (EU)

US PAT 6,013,096; US PAT 6,746,473 For additional US and International patents and patent pending information go to [www.erchonia.com](http://www.erchonia.com).



**ERCHONIA**  
World Leaders in Low Level Laser Technology