

ZERONA®-Z6

ADJUSTABLE
HEIGHT

EASY TO USE TOUCH
SCREEN CONTROL

FLEX ARMS CONFORM
TO VIRTUALLY ANY
POSITION



CARBON
FIBER
FINISH



ZERONA®-Z6 by Erchonia®

The Erchonia Zerona-Z6 laser is specially designed for non-invasive body contouring. The laser is applied externally and has been proven through double-blind, multi-site, placebo controlled clinical trials to emulsify adipose tissue.

As with all great technology, there are advancements and the Zerona Z-6 is just that. Through customer feedback and our endless quest to provide the best non-invasive procedure possible, we have now upgraded the Zerona technology to maximize coverage and offer a better overall treatment. With the reconfiguration of the diodes on the new Z6, you can now be assured your patients are receiving the most researched and clinically effective body contouring procedure available.

The Zerona-Z6 Advantage:

- ▶ **No pain, no downtime and no risk**
- ▶ **Proven effective**
- ▶ **Superior treatment area coverage**



SPECIFICATIONS

Configuration: 6 Line Generated Class 2 Laser Diode Modules

Wavelength: 635nm

Modulation: Constant Wave (CW)

Display: Full Color TFT Touch Screen Control Center

Adjustments:

- 1) 44" (111.76 cm) Vertical Arm Height Adjustment.
- 2) Four Independent Adjustable Arms for Desired Laser Concentration

Power Source: 100-240VAC 50-60Hz

Chassis:

- 1) Metal Frame Powder Coated for Ease of Cleaning
- 2) 4 Locking Anti-Static Casters

Housing: Black Carbon Fiber Finish Thermoformed from Non-Allergen Material

Weight: 71.lbs / 32. kg

Accessories: 2-Keys, Laser Safety Glasses

Compliant to Quality / ISO 13485 - Medical Device / IEC60825-1 Laser Safety, IEC 60601-1 Safety, IEC 60601-1-2 EMC, CE Mark, CB Cert

FDA Laser Class 2 / FDA Device Class II / EU Device Class IIa, Laser Class II

Two Year Full Warranty

US PAT 6,013,096; US PAT 6,746,473; US PAT 8,439,927; US PAT 8,366,756; US PAT 8,348, 985; US PAT 8,097,029; US PAT 7,922,751; US PAT 8,932,338. For additional US and International patents and patent pending information go to www.erschonia.com.

www.myzerona.com
877-699-3766 Toll Free

