

## Introducing the Newest Line of Cold Lasers from Erchonia®

The EVRL Violet Laser System is the newest Erchonia handheld laser in "smart phone" size with a user-friendly graphical user interface (GUI) and a sleek new cord free design. Operating at a wavelength of 405nm and 640nm, the EVRL Laser provides the practitioner with an innovative tool for treatment versatility.

## The EVRL Laser System Advantage:

- · Anti-bacterial & Anti-microbial
- Red & Violet Laser Design
- Touchscreen GUI
- · Made in the USA



## **SPECIFICATIONS**

Configuration: (1) <10 mW Line Generated Red Laser Diode, (1) <5 mW Line Generated Violet Laser Diode with Patented Optics

Wavelengths: 640nm & 405nm

Modulation: Constant Wave (CW) & Pulsed (Hz) (Red Only)

40 User Defined Memory Channels

Display: Full Color TFT Touchscreen Technology for Ease of Use

Power Source: 3.7 VDC Lithium-ion Polymer Battery

Enclosure: Micro-Textured Anodized Aircraft Aluminum for Strength and Ease of Cleaning

Weight: 0.95lbs (.42 kg) (Device Only)

Accessories: Power Supply, Wireless Charging Base, Carrying Case and Belt Clip

Compliant to ISO 13485 Medical Device Quality, IEC 60601-1 Safety, 60601-1-2 EMC, IEC 60825-1 Laser Safety, CE Mark

888.242.0571 Toll Free

www.erchonia.com

