

FX 635



INTRODUCING ERCHONIA'S FX 635 LASER

The world's only laser designed and FDA Market Cleared to relieve chronic low back pain of the musculoskeletal origin and relieve the pain and inflammation of the heel associated with plantar fasciitis. Our FDA market cleared laser has successfully been proven in multiple Level (1) clinical trials to effectively treat these painful conditions. The pain-relieving, healing powers of this low-level laser technology can reduce inflammation and eliminate pain to help your patients return to their active lifestyle in only a few weeks time.

Unlike other lasers on the market, the FX 635 is proven to promote cellular function through painless bio-stimulation. The revolutionary design of the FX 635 Laser allows for dependable positioning of the three patented laser diode arms perfectly to target these pain conditions. It provides patients experiencing chronic pain with a non-invasive, efficacious solution to quickly improve comfort and quality of life without surgery or harmful NSAIDs or opioids. While other lasers may claim to offer relief, only the FX 635 has been FDA market cleared for these conditions.

SPECIFICATIONS:

Configuration: (3) Class 2 17.25mW Line Generating Diodes

Wavelength: 635nm

Modulation: Variable

Weight: 65lbs (29.48kg)

Height: 70in (177.8 cm) (Average - Adjustable)

Two Independent Adjustable Arms For Desired Laser Concentration

4 Locking Anti-Static Casters

Display: Full Color Touch Screen Control

Power Source: 100-240VAC 50-60Hz

Chassis: Powder Coated Aircraft Aluminum

Housing & Covers: Non-Allergenic Material, Spray Finished

Accessories: 2-Keys, Power Cord, Laser Safety Glasses

Compliant to: ISO 13485 Medical Device Quality, IEC 60825-1 Laser Safety, IEC

60601-1-2 EMC, IEC 60601-1 Safety, CE Mark, CB Certified

FDA Laser Class 2 / FDA Device Class II / EU Device Class IIa, Laser Class II



www.erchonia.com
(866) 551-7847 Toll Free

