

FDA Cleared - To Target Onychomycosis



BY ERCHONIA
LunulaLaser™

Provided by: Erchonia Corporation.
For information call: (855) 458-6852 (4LUNULA) or visit: www.Lunula.com



Healthy, Clear Nails Made Simple

The 1st and only non-thermal, touchless, painless treatment of Onychomycosis to target the fungus and stimulate new clear nail growth.

SPECIFICATIONS

Configuration: (2) Line Generated Class 2 Laser Diode Modules
Wavelength: 635nm & 405nm
Modulation: Constant Wave (CW)
Display: Full Color TFT Touch Screen Control Center
Power Source: 100-240VAC 50-60Hz
Chassis: Powder Coated Aircraft Aluminum
Treatment Area: 304 Stainless Steel for Anti-Microbial Properties
Housing: Injection Molded Process with Non-Allergen Material/Plastic
Weight: 23lbs. (10.43 kg)
Accessories: 2-Keys, Power Cord & Laser Safety Glasses
Compliant to: ISO 13485 Medical Device Quality; IEC 60825-1 Laser Safety; IEC 60601-1-2 EMC; IEC 60601-1 MEE Safety; CE Mark; CB Cert; CMDCAS
Laser Class 2: Device Class: II (USA); 2a (EU)



*Individual Results May Vary.

US PAT 6,013,096; US PAT 6,746,473; US PAT 8,409,264; US PAT 8,814,924; US PAT 8,097,029;
US PAT 7,118,588; US PAT 7,947,067; US PAT 7,922,751 and several U.S. and International Patents Pending.

Lunula Laser DataSheet_Rev 3 05/17



Introducing the Erchonia® Lunula Laser®, the revolutionary low level laser therapy bringing new hope to people suffering from onychomycosis. Safe and effective, Lunula Laser is the first and only non-thermal laser to receive FDA 510(k) Market Clearance for Onychomycosis. Lunula Laser poses none of the risks and harmful side effects of oral anti-fungal medications and is painless, unlike other laser therapies.

The Erchonia® Lunula Laser®, uses low level laser light to target onychomycosis and requires very little time or set-up for physicians or their staff. It's the only system you can turn on and walk away—the device already has a pre-set protocol. In fact, as little as (4) 12-minute sessions are needed to treat the affected area.*

The Lunula Advantage:

- No pain, no downtime and no risk
- Avg. 6.1 clearance at 6 months*
- No dangerous smoke plumes
- Unattended procedure
- 89% of patients respond
- Made in the USA



To learn more call or visit us online today!
(855) 458-6852 (4LUNULA) www.lunula.com