

VLS

Powerful Handheld Technology

Easy to Use
Touchscreen
Technology



Introducing the Newest Line of Cold Lasers from Erchonia®

The Vet Laser System is a revolutionary new hand held laser device with a user-friendly graphical user interface (GUI) and a sleek new cord free design. Now you can treat an array of conditions with easy to use preset protocols, or user programmable channels.

SPECIFICATIONS

Configuration: (2) 7.5mW Line Generated Laser Diodes with Patented Optics

Wavelength: 635nm

Modulation: Constant Wave (CW) & Pulsed (Hz)

Display: Full Color TFT Touchscreen Technology for Ease of Use

Power Source: 3.7 VDC Lithium-ion Polymer Battery

Enclosure: Micro-Textured Anodized Aircraft Aluminum for Strength and Ease of Cleaning

Weight: Less Than 1lb. (.45 kg) (Device Only)

Accessories: Power Cord, Wireless Charging Base, Carrying Case and Belt Clip

Compliant to ISO 13485 Medical Device Quality, EN/ISO/IEC 60601-1, 60601-1-2, IEC 60825-1, CE Mark, CMDCAS

Laser Class 2 / Device Class II (USA); 2a (EU)

US PAT 6,013,096; US PAT 6,746,473. For Additional US and International Patents and Patent Pending Information go to www.erschonia.com



ERCHONIA®

World Leader in Low Level Laser Technology

888.242.0571 Toll Free

www.laserpetcare.com